

## Rodent Shavers

In many experimental applications it is recommended that the animal is shaved prior to surgery. Since rodents have very thin hair not all electrical shavers can be used successfully, but these two shavers have proven to work very well with this application.

### ISIS

The Isis shaver is a very small high quality and flexible shaver delivered with a 2,4cm blade. With this blade the shaver can cut as close as 0,1mm to the skin. Thanks to the small footprint the Isis can be located very close to where the surgery takes place.

The Isis shaver can only be used without power cord, and stands firmly in its charger stand while charging the battery. An indicator light is illuminated while charging. The blade can be replaced when needed.

- Ergonomic
- Rechargeable
- Very quiet and smooth
- Blade only 2,4cm wide
- Weight only 110 grams
- Stand-up charger with indicator light
- Operates only without cord

### ANDIS

The Andis shaver is a small and shaver delivered with a 2,7cm blade. With this blade the shaver can cut as close as 0,1mm to the skin. Due to the small footprint also the Andis can be located very close to where the surgery takes place.

The Andis shaver can only be used



Left: Isis Right: Andis

without power cord, and stands firmly in its charger stand while charging the battery. An indicator light is illuminated while charging. The blade can be replaced when needed.

- Rechargeable
- Blade 2,7cm wide
- Weight 140 grams
- Stand-up charger with indicator light
- Operates only without cord

### Ordering Info

Art.No. GT420 Isis rodent shaver  
Art.No. 3515 Andis Shaver, 2,7cm blade